


공 컬럼

XK Column



	제품명	대응 컬럼	포장	코드 번호	재질	그림 속 번호
패킹 관련 제품	Packing reservoir RK 16/26 (XK 16 및 XK 26 용)	XK16, K26	1개	28989858		
	Packing reservoir RK 50 (XK 50용)	XK50	1개	28989861		
	Packing connector XK 16	XK16	1개	18115344		
	Packing Connector XK 26	XK26	1개	18115345		
	Bottom piece XK 16	XK16	1개	28989885	A	1
바텀 피스	Bottom piece XK 26	XK26	1개	28989887	A	1
	Bottom piece XK 50	XK50	1개	28989888	A	1
O-링	O-ring (내경16 mm컬럼용)	XK16	5개	19016301	E	2
	O-ring (내경26 mm컬럼용)	XK26	5개	28978227	E	2
	O-ring (내경50 mm컬럼용)	XK50	5개	28978228	E	2
플런저	Plunger (XK 16, AC 16용)	XK16	1개	18103180	B	3
	Plunger (XK 26, AC 26용)	XK26	1개	18103181	B	3
	Plunger (XK 50용)	XK50	1개	18875801	B	3
서포트 스크린	Support screen (내경16 mm컬럼용)	XK16	5개	19065101	D	4
	Support screen (내경26 mm컬럼용)	XK26	5개	18937701	D	4
	Support screen (내경50 mm컬럼용)	XK50	5개	19066401	D	4
넷 링 10 µm	Net ring 10 µm (내경16 mm컬럼용)	XK16	5개	18876101	B,C	5
	Net ring 10 µm (내경26 mm컬럼용)	XK26	5개	18876001	B,C	5
	Net ring 10 µm (내경50 mm컬럼용)	XK50	5개	18875901	B,C	5
글래스 튜브	Glass tube 20 cm (XK 16/20용)	XK16/20	1개	19031501	F	6
	Glass tube 40 cm (XK 16/40용)	XK16/40	1개	19011301	F	6
	Glass tube 70 cm (XK 16/70용)	XK16/70	1개	19011401	F	6
	Glass tube 100 cm (XK 16/100용)	XK16/100	1개	19011501	F	6
	Glass tube 20 cm (XK 26/20용)	XK26/20	1개	18100084	F	6
	Glass tube 40 cm (XK 26/40용)	XK26/40	1개	19014501	F	6
	Glass tube 70 cm (XK 26/70용)	XK26/70	1개	19014601	F	6
	Glass tube 100 cm (XK 26/100용)	XK26/100	1개	19014701	F	6
	Glass tube 20 cm (XK 50/20용)	XK50/20	1개	18100085	F	6
	Glass tube 30 cm (XK 50/30용)	XK50/30	1개	19132601	F	6
	Glass tube 60 cm (XK 50/60용)	XK50/60	1개	19052501	F	6
	Glass tube 100 cm (XK 50/100용)	XK50/100	1개	19050901	F	6
자켓	Acrylic jacket 20 cm (XK 16/20용)	XK16/20	1개	18001172	G	7
	Acrylic jacket 40 cm (XK 16/40용)	XK16/40	1개	18001173	G	7
	Acrylic jacket 70 cm (XK 16/70용)	XK16/70	1개	18001174	G	7
	Acrylic jacket 100 cm (XK 16/100용)	XK16/100	1개	18001175	G	7
	Acrylic jacket 20 cm (XK 26/20용)	XK26/20	1개	18100082	G	7
	Acrylic jacket 40 cm (XK 26/40용)	XK26/40	1개	18001167	G	7
	Acrylic jacket 70 cm (XK 26/70용)	XK26/70	1개	18001168	G	7
	Acrylic jacket 100 cm (XK 26/100용)	XK26/100	1개	18879101	G	7
	Acrylic jacket 20 cm (XK 50/20용)	XK50/20	1개	18100083	G	7
	Acrylic jacket 30 cm (XK 50/30용)	XK50/30	1개	18001178	G	7
	Acrylic jacket 60 cm (XK 50/60용)	XK50/60	1개	18001177	G	7
	Acrylic jacket 100 cm (XK 50/100용)	XK50/100	1개	18876501	G	7
어댑터	Adapter XK 16	XK16	1개	28989876		8
	Adapter XK 26	XK26	1개	28989877		8
	Adapter XK 50	XK50	1개	28989880		8
엔드 피스	End piece XK 16	XK16	1개	18648801		9
	End piece XK 26	XK26	1개	18648901		9
	End piece XK 50	XK50	1개	18879701		9
기타	Stop plug, 1/16" female		5개	11000464		

주) 모든 XK 컬럼에는 AK 어댑터가 1 개씩 장착되어 있습니다.

재 질	A Reinforced acetal plastic	E Ethylenepropylenedione polymer (EPDM)
	B Polypropylene(reinforced glassfibre)	F Borosilicate glass
	C Polyamide(Nylon)	G Acrylic plastic
	D Polypropylene	H Polypropylene with glass fibre reinforcement (PP-GF)